

*Innovation  
with  
Diversity*

# Aspirator

**D-30**

**D-50**

**D-90**

Emergency Dept.

Operation Theater

General Surgery

Neurosurgery

Liposuction

Intensive Care Unit

Gynaecology



CE 0123

## D-90

### TECHNICAL FEATURES

<ul style="list-style-type: none"><li>▶ Motor</li><li>▶ Power Feeding</li><li>▶ ISO 10091-1 Classification</li><li>▶ Max Vacuum (adjustable)</li><li>▶ Max free air flow rate</li><li>▶ Noise Level</li><li>▶ Power consumption</li><li>▶ Fuse</li><li>▶ Duty cycle</li><li>▶ Foot Switch</li><li>▶ Outside Dimension</li><li>▶ Year of Warranty</li><li>▶ Collection Jar</li><li>▶ Optional</li></ul>	<p>Reciprocating Oil less piston pump 220V-50Hz, 110V-50Hz HIGH VACUUM -0.90 Bar, -90kpa, -675 mmHg 90 L/min ≤45dB(A) 250VA 1 x 2.0 A 250V Non-stop operation, Explosion proof pedal switch 970 x 565 x 475 mm 1 2 x 4000ml or 2 x 2000ml (Polycarbonate / Polysulfone) Autoclavable Triple Flow Safety, 10x Anti-Bacterial Filter, Vacuum Silicon Cup, 5000ml Jar, Disposable Liquid Bag, Change Over Switch, ≤40 dB(A)</p>
--	--

## D-50

### TECHNICAL FEATURES

<ul style="list-style-type: none"><li>▶ Motor</li><li>▶ Power Feeding</li><li>▶ ISO 10091-1 Classification</li><li>▶ Max Vacuum (adjustable)</li><li>▶ Max free air flow rate</li><li>▶ Noise Level</li><li>▶ Power consumption</li><li>▶ Fuse</li><li>▶ Duty cycle</li><li>▶ Foot Switch</li><li>▶ Outside Dimension</li><li>▶ Year of Warranty</li><li>▶ Collection Jar</li><li>▶ Optional</li></ul>	<p>Reciprocating Oil less piston pump 220V-50Hz, 110V-50Hz HIGH VACUUM -0.90 Bar, -90kpa, -675 mmHg 50 L/min ≤45dB(A) 250VA 1 x 2.0 A 250V Non-stop operation Explosion proof pedal switch 970 x 565 x 475 mm 1 2 x 4000ml or 2 x 2000ml (Polycarbonate / Polysulfone) Autoclavable Triple Flow Safety, 10x Anti-Bacterial Filter, Vacuum Silicon Cup, 5000ml Jar, Disposable Liquid Bag, Change Over Switch, ≤40 dB(A)</p>
--	---

## D-30

### TECHNICAL FEATURES

<ul style="list-style-type: none"><li>▶ Motor</li><li>▶ Power Feeding</li><li>▶ ISO 10091-1 Classification</li><li>▶ Max Vacuum (adjustable)</li><li>▶ Max free air flow rate</li><li>▶ Noise Level</li><li>▶ Power consumption</li><li>▶ Fuse</li><li>▶ Duty cycle</li><li>▶ Foot Switch</li><li>▶ Outside Dimension</li><li>▶ Year of Warranty</li><li>▶ Place of Manufacturing</li><li>▶ Collection Jar</li><li>▶ Optional</li></ul>	<p>Reciprocating Oil less piston pump 220V-50Hz, 110V-50Hz HIGH VACUUM -0.85 Bar, -85kpa, -650 mmHg 30 L/min ≤45dB(A) 250VA 1 x 2.0 A 250V Non-stop operation Explosion proof pedal switch 970 x 500 x 465 mm 1 2 x 1000ml or 2 x 2000ml (Polycarbonate / Polysulfone) Autoclavable Triple Flow Safety, 10x Anti-Bacterial Filter, Vacuum Silicon Cup, 3000ml Jar, Disposable Liquid Bag, Change Over Switch, ≤40 dB(A)</p>
---	---

DISTRIBUTOR:

A Member of:



### CONTACT:

DISPOMEDICAL UK LTD  
21 Duckworth Lane,  
Bradford BD95ER UK

Watsapp : +962776663777

Phone : +441133452242

Email : [export@dispomedical.co.uk](mailto:export@dispomedical.co.uk)

Website : [www.dispomedical.co.uk](http://www.dispomedical.co.uk)