




**D-1000 V**  
Best choice for veterinary

### Outstanding Clinical Application:

D-1000V has excellent performance and different applications according to user's demands Felidae, Canine, Ovine, Equine, Bovine/buffaloes, Swine, User-defined.

### Applications:

Abdominal, Vascular, Reproduction, Musculoskeletal, Cardiac



### Advanced imaging technologies:

- Exclusive employ double wave form imaging technology
- Speckle noise reduction technology
- Tissue specific imaging technology
- Tissue harmonic imaging technology
- Pulse wave imaging mode (PW)
- Trans abdominal application
- Phase Inversion Harmonic Imaging
- One -click automatic image optimization



### Specific Veterinary Transducers



Convex

Linear

Micro-convex

Vet Micro-convex

Linear Rectal

### DISPOMEDICAL UK LTD

21 Duckworth Lane, Bradford BD9 5ER UK

Phone: +44(0)8450779337 Email: info@dispomedical.co.uk

## Features:

- Outstanding portable model compact and streamlined design
- 12.1- inch high resolution medical monitor with 90 degrees rotation
- Two activated transducer connectors
- Auto-freezing function
- Imaging function: B zoom
- One key full screen function
- Abundant measurement software packages
- Auto generation BMP format report
- Support USB import and export
- DICOM 3.0
- Review: Cine-loops, 1024 frames
- Built in lithium battery
- Optional trolley and suitcase

## Clinic Value and Applications




Canine Liver




Canine Intestine



Canine Kidney



Bovine Gender



Bovine Follicle



Canine Pregnancy

## DISPOMEDICAL UK LTD

21 Duckworth Lane, Bradford BD9 5ER UK

Phone: +44(0)8450779337 Email: info@dispomedical.co.uk



**DISPOMEDICAL UK LTD**

21 Duckworth Lane, Bradford BD9 5ER UK

Phone: +44(0)8450779337 Email: [info@dispomedical.co.uk](mailto:info@dispomedical.co.uk)